


KN95 PROTECTIVE MASK

Breathable filtration, layers of protection





- **Item Name :** KN95
- **Product structure :** It is made up of mask body ,nose clip and mask belt
- **Packing :** 10PCS/Box
- **Production Company :** Shantou Guoheng Industrial Co.,Ltd
- **Address :** Chendian Industrial Zone , Chaonan District , Shantou , Guangdong , China
- **Made in China**





Instructions:

- **Product use and scope of application :** Applicable to the protection of non-oil particles , effective protection against PM2.5、 dust and haze , blocking droplets and filtering the air in a variety of microorganisms
- **Notes, warning and reminders :**
 - ① This product can not be used in anoxic environment or in patients with breathing difficulties
 - ② This product can not be used as a gas mask in toxic gas environment
 - ③ Do not use if the package is damaged or expired
- **Special storage conditions and methods :** Store in a room with relative humidity less than 80% , no corrosive gas and good ventilation to avoid high temperatures
- **Period of validity :** Stored under specified conditions and valid for 24 months



- ☑ Prevent PM2.5
- ☑ Prevent Dust
- ☑ Prevent Bacteria
- ☑ Prevent Haze



KN95 mask wearing method



① Place the nose clip over the mask and hold it against your chin



② After straightening the ear strap, adjust the ear to feel as comfortable as possible



③ Move from the middle of nose clip to both sides and press until the nose clip to fit the bridge of the nose



④ Check the air-leaking of your face and the mask



Sealing Bag (1PCS/Bag)

Size : 115 X 180 mm Material : 3C HDPE Print color : Black



Color box packaging (10PCS/Box)

Size : 170 X 115 X 50 mm Material : 400g white cardboard Print color: Black+gray



Outer Carton (100 boxes / carton)

Size : 590 X 520 X 355 mm Material: Double Wall strengthening Print color: Black

Certificate of Compliance



No. 0H200402S.SGI0073

Technical Construction File no. XMT0202001372LY/PPE

Certificate's
Holder:

Shantou Guoheng Industrial Co., Ltd
Chendian Industrial Zone, Chaonan District,
Shantou, Guangdong, China

Certification ECM
Mark:



Product:

KN95

Model(s):

KN95 (FFP2)

Verification to:

Standard:

EN 149:2001+A1:2009

related to CE Directive(s):

R 2016/425 (Personal Protective Equipment)

Remark: This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the technical documentation received from the manufacturer is satisfactory for the requirements of the ECM Certification Mark. The conformity mark above can be affixed on the products accordingly to the ECM regulation about its release and its use.

Additional information and clarification about the Marking:



The manufacturer is responsible for the CE Marking process, and if necessary, must refer to a Notified Body. This document has been issued on the basis of the regulation on ECM Voluntary Mark for the certification of products. RG01_ECM rev.3 available at: www.entecerma.it

Issuance date: 02 April 2020

Expiry date: 01 April 2025

Reviewer
Technical expert
Amanda Payne



Approver
ECM Service Director
Luca Bedonni




Ente Certificazione Macchine Srl

Via Ca' Bella, 243 - Loc. Castello di Serravalle - 40053 Valsamoggia (BO) - ITALY
☎ +39 051 6705141 📠 +39 051 6705156 ✉ info@entecerma.it 🌐 www.entecerma.it

Test report 04

The products have been certified by authoritative institutions

 180011112242 (2018) 国认监认字(204)号 检验报告 TEST REPORT  Scan is following us Scan checks for authenticity  W202007701E 浙江省轻工业产品质量检验研究院 Zhejiang Light Industrial Products Inspection and Research Institute 国家纺织服装产品质量监督检验中心(浙江) National Textiles and Garment Quality Supervision Inspection Center(Zhejiang)		Zhejiang Light Industrial Products Inspection and Research Institute National Textiles and Garment Quality Supervision Inspection Center(Zhejiang) Test report Number: W202007701E page1/3 <table border="1"> <tr> <td>Name of Customer</td> <td>Shantou Guobeng Industrial Co., LTD</td> <td>Address</td> <td>Chendian Industrial Zone, Chaonan District, Shantou, Guangdong, China</td> </tr> <tr> <td>Manufacturer</td> <td>---</td> <td>Address</td> <td>---</td> </tr> <tr> <td>Sample Information</td> <td colspan="3"> Name of sample: KN95 MASK Characteristics of sample: WHITE Trademark of sample: --- Specification/model: --- Level: KN95 Category of safety specification: --- Art. No.: --- </td> </tr> <tr> <td colspan="4">Above-mentioned information by Customer-supplied</td> </tr> <tr> <td>The sent way of sample</td> <td>Courier</td> <td>Sample quantity</td> <td>30 pieces</td> </tr> <tr> <td>Receiving Date of Sample</td> <td>2020/04/18</td> <td>Test Category</td> <td>Entrusted inspection</td> </tr> <tr> <td>Date of Testing</td> <td colspan="3">2020-04-18~ 2020-04-19</td> </tr> <tr> <td>Rating Requirements</td> <td colspan="3">GB 2626-2006</td> </tr> <tr> <td colspan="4"> Test Summary: See the attached page for the results. </td> </tr> <tr> <td colspan="4">  Date of Approval: 2020.4.20 </td> </tr> <tr> <td>Remarks</td> <td colspan="3">The sample is not pretreated and the result is for reference.</td> </tr> <tr> <td colspan="4"> Approved by: 曹丽勤 </td> </tr> </table>		Name of Customer	Shantou Guobeng Industrial Co., LTD	Address	Chendian Industrial Zone, Chaonan District, Shantou, Guangdong, China	Manufacturer	---	Address	---	Sample Information	Name of sample: KN95 MASK Characteristics of sample: WHITE Trademark of sample: --- Specification/model: --- Level: KN95 Category of safety specification: --- Art. No.: ---			Above-mentioned information by Customer-supplied				The sent way of sample	Courier	Sample quantity	30 pieces	Receiving Date of Sample	2020/04/18	Test Category	Entrusted inspection	Date of Testing	2020-04-18~ 2020-04-19			Rating Requirements	GB 2626-2006			Test Summary: See the attached page for the results.				 Date of Approval: 2020.4.20				Remarks	The sample is not pretreated and the result is for reference.			Approved by: 曹丽勤			
Name of Customer	Shantou Guobeng Industrial Co., LTD	Address	Chendian Industrial Zone, Chaonan District, Shantou, Guangdong, China																																																
Manufacturer	---	Address	---																																																
Sample Information	Name of sample: KN95 MASK Characteristics of sample: WHITE Trademark of sample: --- Specification/model: --- Level: KN95 Category of safety specification: --- Art. No.: ---																																																		
Above-mentioned information by Customer-supplied																																																			
The sent way of sample	Courier	Sample quantity	30 pieces																																																
Receiving Date of Sample	2020/04/18	Test Category	Entrusted inspection																																																
Date of Testing	2020-04-18~ 2020-04-19																																																		
Rating Requirements	GB 2626-2006																																																		
Test Summary: See the attached page for the results.																																																			
 Date of Approval: 2020.4.20																																																			
Remarks	The sample is not pretreated and the result is for reference.																																																		
Approved by: 曹丽勤																																																			
Number: W202007701E page2/3 <table border="1"> <thead> <tr> <th>ITEM</th> <th>STANDARD (KN95)</th> <th>RESULT</th> <th>RATING</th> </tr> </thead> <tbody> <tr> <td colspan="4">1. FILTER EFFICIENCY(SALT MEDIUM) (GB2626-2006)(%)</td> </tr> <tr> <td>-</td> <td>≥95.0</td> <td>99.3</td> <td>PASS</td> </tr> <tr> <td colspan="4">2. RESPIRATORY RESISTANCE (GB2626-2006)(Pa)</td> </tr> <tr> <td>-INSPIRATORY RESISTANCE</td> <td>≤350</td> <td>61</td> <td>PASS</td> </tr> <tr> <td>-EXPIRATORY RESISTANCE</td> <td>≤250</td> <td>56</td> <td>PASS</td> </tr> </tbody> </table> 		ITEM	STANDARD (KN95)	RESULT	RATING	1. FILTER EFFICIENCY(SALT MEDIUM) (GB2626-2006)(%)				-	≥95.0	99.3	PASS	2. RESPIRATORY RESISTANCE (GB2626-2006)(Pa)				-INSPIRATORY RESISTANCE	≤350	61	PASS	-EXPIRATORY RESISTANCE	≤250	56	PASS	Number: W202007701E page3/3 Test report Picture(s) of sample  —End of report— 																									
ITEM	STANDARD (KN95)	RESULT	RATING																																																
1. FILTER EFFICIENCY(SALT MEDIUM) (GB2626-2006)(%)																																																			
-	≥95.0	99.3	PASS																																																
2. RESPIRATORY RESISTANCE (GB2626-2006)(Pa)																																																			
-INSPIRATORY RESISTANCE	≤350	61	PASS																																																
-EXPIRATORY RESISTANCE	≤250	56	PASS																																																

Test result

DISPOSABLE MEDICAL MASK

一次性使用医用口罩

Breathable filtration, layers of protection



Item Name: Disposable Medical Mask

Size: 17.5 X 9.5 cm

Packing: 50PCS /Box

Material: Double-layer non-woven with melt-blown non-woven filter layer

Carried Standard: YY/T 0969-2013

Production license number: Guangdong-shantou medical device production emergency

Record 20200013


Medical device registration certificate number: Guangdong-shantou medical device emergency

Record 20200023

Production Company : Shantou Guoheng Industiral Co.,Ltd

Address: Chendian Industrial Zone , Chaonan District , Shantou , Guangdong , China





Risk of Infection

Masks can filter out certain pollutants, but do not eliminate the risk of infection



Timely replacement

When the mask is broken or contaminated, please replace it in time



Sensitive skin

Sensitive skin may be allergic to the material, it is recommended to test the material before use



Risk of Infection

It might spread disease if you use it in wrong way




Fit degree

Please check the air tightness every time you wear it to ensure the mask fits the face



Pregnant women

Pregnant women and people with respiratory problems should use it on the advice of a doctor





Effective against

Bacteria

Droplet spread

Peculiar smell

Dust



Right way



Flatten the mask and hang it on both sides of your ears



Fix the nose clip in the shape of the nose to prevent dirty air from entering



Please pull the mask to your chin to present stereoscopic breathing space



Sealing Bag (10PCS/Bag)

Size : 120 X 230 mm Material : 3C HDPE Print color : Blue+Black



Color box packaging (50PCS/box)

Size : 190 X 85 X 105 mm Material : 400g white cardboard Print color: Blue+gray



Outer Carton (40 boxes / carton)

Size : 550 X 410 X 360 mm Material: Double Wall strengthening Print color : Black



The 3rd Party Certificate of FDA Medical Device Registration & Listing

Note: This doc is Not issued by FDA. We XMT create it to make easier for customer to present information.
The following registered establishment & device information, are excerpted from FDA database.

Establishment

SHANTOU GUOHENG INDUSTRIAL CO., LTD

No. 5, North TongFu Road, ShuiKou, HuiCheng,
Shantou Guangdong, CN 515152

Registration Number: Awaiting Assignment | FEI Number: Awaiting Assignment

Status: **Active** | Date of Registration Status: **2020**

{ More details [*a](#) }

Owner/Operator

Shantou Guoheng Industrial Co., Ltd

No. 5, North TongFu Road, ShuiKou, HuiCheng,
Shantou Guangdong, CN 515152

Owner/Operator Number: **10063893**

{ More details [*b](#) }

Official Correspondent

Weijian Chen

Phone: 086-754-84416886

US Agent

Jiali Chen

2010 3rd Ave, #2, Los Angeles, CA US 90018

Phone: (+1) 415 2977106 | E-mail: jchen01234@gmail.com

Devices Listing Information

#01	Listing Number:	D37****	{ More details *01 }
	Proprietary Name:	Disposable Medical Mask	
#02	Listing Number:	N/A	{ More details *02 }
	Proprietary Name:	N/A	
#03	Listing Number:	N/A	{ More details *03 }
	Proprietary Name:	N/A	
#04	Listing Number:	N/A	{ More details *04 }
	Proprietary Name:	N/A	
#05	Listing Number:	N/A	{ More details *04 }
	Proprietary Name:	N/A	

Caution: Please careful protect your device Listing Number.

* Remark: More devices continued on next page [if applicable].

--- Above lists one (1) devices of this establishment has been listed with U.S. FDA ---



Professional FDA register services, from XMTEST www.xmtest.org
Need help? Contact us, at +86 021-6104 9808 or xmtest@xmtest.org

Reg Agent

CE Documentation Review



No. 0H200325M.SGI0243

Holder: Shantou Guoheng Industrial Co., Ltd.
Chendian Industrial Zone, Chaonan District,
Shantou, Guangdong, China

Review goal: Verification of the presence of Technical Documentation compatible with the Medical Devices Directive 93/42/EEC Annex VII

Product: Disposable Medical Mask (Not Sterile)

Model(s): GY-L, GY-S

Classification: Class I (Not Sterile)
(accordingly to the Manufacturer's declaration)

Review output: This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the Technical Documentation shared with us by the manufacturer is compatible with the European Standard for Medical Devices.
The manufacturer is responsible for the CE Marking process, and not exempted to carry out all necessary compliance activities. This document has been issued on the basis of the regulation on ECM Voluntary Mark for the certification of products. RG01_ECM rev.3 available at: www.entecerma.it

Date of issue 25 March 2020

Expiry date 24 March 2025

Approver
ECM Service Director
Luca Redonni



Technical Expert
Amanda Pavina



Ente Certificazione Macchine

Via Cà Bella, 243 - 40053 Valsamoggia Loc. Castello di Serravalle (Bo) Italy
☎ +39.0516705141 📠 +39.0516705156 ✉ info@entecerma.it 🌐 www.entecerma.it

检验结果

检验结果

The products have been certified by authoritative institutions



深圳市计量质量检测研究院 Shenzhen Academy of Metrology & Quality Inspection MA Metrology A Inspection TEST REPORT No.: WT204013755 Page 1 of 4 Pages Sample Description: Disposable medical mask Model/Specification/Grade: Applicant: Shantou Guoheng Industrial Co.,LTD Applicant Address: Chendian Industrial Zone,Chaonan District,Shantou,Guangdong,China Date of Receipt: 2020-02-25 Test Period: 2020-02-25 to 2020-03-03 Test Location: Longhua Experimental Base Shenzhen Academy of Metrology & Quality Inspection (Stamp) Approved by: 滕万红(副部长) Issue Date: 2020-03-23 Signature: 深圳市计量质量检测研究院 Shenzhen Academy of Metrology & Quality Inspection http://www.smqc.com.cn 电子邮箱(E-mail): kfy@smqc.com.cn CMAE 19197302 & 201719001402 龙华实验基地: 深圳龙华区民治大道民康路114号 咨询电话: 0755-27529955 16位: 0755-27528707 邮编: 518111 Longhua Experimental Base No.114, Wukong North Road, Minzhi Avenue Longhua District, Shenzhen Tel:0755-27529955	深圳市计量质量检测研究院 Shenzhen Academy of Metrology & Quality Inspection TEST REPORT No.: WT204013755 Page 2 of 4 Pages Sample Information: Sample Description: Disposable medical mask Trade Mark: Model/Specification/Grade: Serial/Batch No. of Sample: Manufactured Date: Manufacturer: Shantou Guoheng Industrial Co.,LTD Manufacturer Address: Sample Quantity: 50PCS Sample Description before Testing: Normal. Client Information: Applicant: Shantou Guoheng Industrial Co.,LTD Applicant Address: Chendian Industrial Zone,Chaonan District,Shantou,Guangdong,China Applicant Telephone: Applicant Post Code: Test Information: Date of Receipt: 2020-02-25 Applicant No.: 8227992 Environment Condition: (18~25)℃ (30~70)%RH Sampling Method: Delivered by Applicant Judgment Basis: YY/T 0969-2013 Test Standard: YY/T 0969-2013and other method standards refer to next pages Test Conclusion: Test result refer to next pages. Tested by: 林彬 林彬 Checked by: 王燕归 王燕归																																				
深圳市计量质量检测研究院 Shenzhen Academy of Metrology & Quality Inspection TEST REPORT No.: WT204013755 Page 3 of 4 Pages <table><tr><th>Test Item</th><th>Requirement</th><th>Test Result</th><th>Conclusion</th></tr><tr><td>1.Bacterial filtration efficiency (BFE)(%) (YY/T 0969-2013)</td><td>≥95</td><td>(YY 0469-2011) 1# 2# 3# 100 100 100</td><td>Conformity</td></tr><tr><td>2.Microorganisms (YY/T 0969-2013)</td><td></td><td>(GB 15979-2002)</td><td>Conformity</td></tr><tr><td>Total amount of bacterial colony (CFU/g)</td><td>≤100</td><td><4</td><td></td></tr><tr><td>Coliform group</td><td>No detected</td><td>Not detected</td><td></td></tr><tr><td>Pseudomonas aeruginosa</td><td>No detected</td><td>Not detected</td><td></td></tr><tr><td>Staphylococcus aureus</td><td>No detected</td><td>Not detected</td><td></td></tr><tr><td>Streptococcus hemolyticus</td><td>No detected</td><td>Not detected</td><td></td></tr><tr><td>Fungal colony</td><td>No detected</td><td>Not detected</td><td></td></tr></table> 附注: 1.细菌过滤效率项目分包给深圳市药品检验研究院(深圳市医疗器械检测中心) (资质认定证书编号:201719001159)。	Test Item	Requirement	Test Result	Conclusion	1.Bacterial filtration efficiency (BFE)(%) (YY/T 0969-2013)	≥95	(YY 0469-2011) 1# 2# 3# 100 100 100	Conformity	2.Microorganisms (YY/T 0969-2013)		(GB 15979-2002)	Conformity	Total amount of bacterial colony (CFU/g)	≤100	<4		Coliform group	No detected	Not detected		Pseudomonas aeruginosa	No detected	Not detected		Staphylococcus aureus	No detected	Not detected		Streptococcus hemolyticus	No detected	Not detected		Fungal colony	No detected	Not detected		深圳市计量质量检测研究院 Shenzhen Academy of Metrology & Quality Inspection TEST REPORT No.: WT204013755 Page 4 of 4 Pages 2.Photo(s) of the sample(s).  END OF REPORT
Test Item	Requirement	Test Result	Conclusion																																		
1.Bacterial filtration efficiency (BFE)(%) (YY/T 0969-2013)	≥95	(YY 0469-2011) 1# 2# 3# 100 100 100	Conformity																																		
2.Microorganisms (YY/T 0969-2013)		(GB 15979-2002)	Conformity																																		
Total amount of bacterial colony (CFU/g)	≤100	<4																																			
Coliform group	No detected	Not detected																																			
Pseudomonas aeruginosa	No detected	Not detected																																			
Staphylococcus aureus	No detected	Not detected																																			
Streptococcus hemolyticus	No detected	Not detected																																			
Fungal colony	No detected	Not detected																																			

Test result

The products have been certified by authoritative institutions

[illegible]

Test result

Test report 04

The products have been certified by authoritative institutions

**深圳市计量质量检测研究院**
Shenzhen Academy of Metrology & Quality Inspection

TEST REPORT



No. WT204016164

Page 1 of 3 Pages

Sample Description: Disposable Medical Mask

Model/Specification/Grade:

Applicant: Shantou Guoheng Industrial Co.,LTD

Applicant Address: Chendian Industrial Zone, Chaonan District, Shantou, Guangdong, China

Date of Receipt: 2020-03-27

Test Period: 2020-03-27 to 2020-04-01

Test Location: Longhua Experimental Base

**Shenzhen Academy of Metrology & Quality Inspection (Stamp)**

Approved by: 何行月

Issue Date: 2020-04-01

Signature: 何行月

深圳市计量质量检测研究院 Shenzhen Academy of Metrology & Quality Inspection <http://www.smq.com.cn>

电子邮箱: smq@smq.com.cn 网站: www.smq.com.cn 地址: 深圳市宝安区西乡街道铁岗社区铁岗路111号 邮编: 518101

龙华实验基地: 深圳市龙华区民治街道民治大道111号 咨询电话: 0755-27528955 传真: 0755-27528957 邮编: 518101

Longhua Experimental Base No.111, Minzhang North Road, Minzhang Avenue, Longhua District, Shenzhen Tel: 0755-27528955

**深圳市计量质量检测研究院**
Shenzhen Academy of Metrology & Quality Inspection

TEST REPORT

No.: WT204016164

Page 2 of 3 Pages

Sample Information:

Sample Description: Disposable Medical Mask

Trade Mark:

Model/Specification/Grade:

Serial/Batch No. of Sample:

Manufactured Date:

Manufacturer: Shantou Guoheng Industrial Co.,LTD

Manufacturer Address: Chendian Industrial Zone, Chaonan District, Shantou, Guangdong, China

Sample Quantity: 200个

Sample Description before Testing: Normal.

Client Information:

Applicant: Shantou Guoheng Industrial Co.,LTD

Applicant Address: Chendian Industrial Zone, Chaonan District, Shantou, Guangdong, China

Applicant Telephone:

Applicant Post Code:

Test Information:

Date of Receipt: 2020-03-27

Applicant No.: 7700078

Environment Condition: (18~25) °C (30~70) %RH

Sampling Method: Delivered by Applicant

Judgment Basis: EN 14683:2019+AC:2019

Test Standard: ISO 22609:2004 and other method standards refer to next pages

Test Conclusion:

Test result refer to next page.

Tested by: 谢丹菊 谢丹菊

Checked by: 陈开江 陈开江

**深圳市计量质量检测研究院**
Shenzhen Academy of Metrology & Quality Inspection

TEST REPORT

No.: WT204016164

Page 3 of 3 Pages

Test Item	Requirement	Test Result	Conclusion
Differential pressure (Pa/cm ²)	Type IIR <60	Average value: 53.7	PASS
		EN 14683: 2019+AC:2019, Annex C	
		Flow rate: 8L/min	
Splash resistance pressure (kPa)	Type IIR ≥16.0	AQL 4.0, the highest corresponding blood pressure: 21.3	PASS
		(Temperature of test environment: 22.3 °C Relative humidity: 81.3%)	
		ISO 22609:2004	



END OF REPORT

Test result

统一社会信用代码
91440514052480165C

(副本)
(副本号:1-1)



扫描二维码登录“国家企业信用信息公示系统”了解更多登记、备案、许可、监管信息。

名称 汕头市国亨实业有限公司

注册资本 人民币壹仟零捌拾玖万元，

类型 有限责任公司(自然人投资或控股)

成立日期 2012年08月28日

法定代表人 蔡标州

营业期限 长期

经营范围

住所 汕头市潮南区陈店镇溪北工业区11号(国亨楼)

生产、加工、销售：针织品、纺织品、乳罩、松紧带、花边、刺绣工艺品、服装、泳衣、内衣及配饰、口罩、纸尿裤、卫生巾。销售：消毒产品、家用电器、电子元件、电话通信设备、照相器材、塑料制品、建筑材料、日用百货、鞋、帽、工艺美术品、文具用品；第二类医疗器械生产；第二类医疗器械经营；货物、技术进出口（法律、行政法规禁止的项目除外；法律、行政法规限制的项目须取得许可后方可经营）。（依法须经批准的项目，经相关部门批准后方可开展经营活动）

登记机关



生产资质认证
获得粤汕械医疗器械、生产备案凭证

关于申办医用口罩应急注册备案和生产的请示

汕头市工业和信息化局：

我司汕头市国亨实业有限公司，是位于汕头市潮南区陈店镇溪北社区一家口罩及防护服的生产企业（统一社会信用代码91440514052480165B）。生产车间位于主楼第4层，符合10万级的洁净车间，并具备微生物试验能力和相关理化试验能力。因目前新型冠状病毒疫情严峻，为响应省政府全办支持疫情防控所需医疗器械供应保障的要求，现依据《广东省药品监督管理局办公室关于一级响应期间对医用口罩等防控急需用品实施特殊管理的补充通知》（粤药监办许【2020】42号）的规定，拟申请办医用外科口罩和防护服应急注册备案和生产，为国家防疫工作贡献企业的一份微薄力量，企业生产的医用口罩、防护服产品，优先满足各级政府部门采购、收藏。

恳请工信部门考虑本地实际情况，支持我司办理医用外科口罩和防护服应急注册备案和生产的诉求，对该情况予以批复。

专此请示，恳请批复。

汕头市潮南区工业和信
息化局

汕头市国亨实业有限公
司

2020年2月23日

抄报：汕头市潮南区工业和信息化局

汕头市潮南区工业和信息化局

汕头市国亨实业有限公
司

2020年2月23日

汕头市防控应急医疗器械备案凭证

备案号：粤汕械应急备 20200023 号

备案人名称	汕头市国亨实业有限公司
备案人组织机构代码	91440514052480165C
备案人注册地址	汕头市潮南区陈店镇溪北工业区11号（国亨楼）
生产地址	汕头市潮南区陈店镇溪北工业区11号（国亨楼）
代理人	***
代理人注册地址	***
产品名称	一次性使用医用口罩（疫情应急产品）
型号/规格	GW-L、GY-S
产品描述	由三层非织造布、口罩带和鼻夹经热压热合制成的一次性使用非灭菌产品。
预期用途	供临床各类人员在非有创操作过程中佩带，覆盖住使用者的口、鼻及下颌，为防止病原体微生物、颗粒物等的直接通过提供一定的物理屏障。
备注	本备案仅在防控新冠肺炎疫情期间有效。
备案单位和日期	汕头市市场监督管理局 备案日期：2020年2月29日
变更情况	

汕头市防控应急医疗器械生产备案凭证

备案编号：粤汕械生产应急备 20200013 号

企业名称	汕头市国亨实业有限公司			
住 所	汕头市潮南区陈店镇溪北工业区11号（国亨楼）			
生产地址	汕头市潮南区陈店镇溪北工业区11号（国亨楼）			
法定代表人	蔡标州	企业负责人	蔡标州	
生产范围	2017分类目录14-13			
生产产品列表	产品名称	产品备案号	登记日期	备注
	一次性使用医用口罩	粤汕械应急备 20200023 号	2020/2/29	本备案仅在防控新冠肺炎疫情期间有效。

备案部门（公章）

备案日期：2020年2月29日

备案专用章

汕头市防控应急医疗器械备案凭证

汕头市国亨实业有限公司：

根据相关法规要求，对你单位申请的防控应急医疗器械：一次性使用医用口罩（疫情应急产品）予以备案，备案号：粤汕械应急备 20200023 号（本备案仅在防控新冠肺炎疫情期间有效）。

汕头市市场监督管理局

备案专用章

日期：2020年2月29日

Thank You !